



FEATURES

- direct measurement of $H_p(10)$ and $H_p(0.07)$ over the entire energy range
- instant readout
- extremely easy readout process
- passive operation
- not sensitive to EM and RF interference
- operation at high dose rates
- operation at pulsed fields
- PTB approved model available

DIS-1

Direct Ion Storage Dosimeter

The RADOS DIS-1 personal dosimeter is based on an ionisation chamber combined with a modern electronic Direct Ion Storage (DIS) memory cell. The Ion chamber is widely used as a reference detector in radiation detection and is now available in everyday dosimetry applications.

The DIS-1 dosimeter could be described as a passive electronic TLD or Film badge, which can be read instantly and non-destructively without any loss of dose information. This unique feature allows the user of the DIS-1 to read his/her accumulated doses whenever necessary.

The DIS-1 dosimeter has a small, lightweight, rugged and watertight construction, which makes the DIS-1 dosimeter reliable and easy to use.

The radiological range of the DIS-1 covers the entire $H_p(10)$ and $H_p(0.07)$ photon and beta energies without any compromises. The wide dose and energy range, the ability to operate in pulsed fields and the performance at high dose rates make DIS-1 an ideal device for all kinds of radiation dosimetry applications. The DIS-1 allows for the detection of heavy, high-energy ions and its immunity to any external interference is unequalled. There are no deviations in the dose readings even at very high EM or RF fields.



health physics

A Mirion Technologies Division


Featuring:

RADOS


TECHNICAL SPECIFICATIONS:

<p>Physical Characteristics</p>	<ul style="list-style-type: none"> • detector type: three TMDIS (Direct Ion Storage) detectors and two MOSFET detectors • sensitive to gamma, X-ray and beta radiation • insensitive to neutrons (<5 %) • instant readout of ICRU dose equivalents: <ul style="list-style-type: none"> - Hp(10) 1 μSv to 40 Sv (0.1 mrem to 4000 rem) ¹⁾ - Hp(0.07) 10 μSv to 40 Sv (1 mrem to 4000 rem) ¹⁾ • calibration accuracy: <ul style="list-style-type: none"> ±5 % at 1 mSv ¹³⁷Cs Hp(10) ±10 % at 10 mSv ¹³⁷Cs Hp(0.07) • energy response in the dose range up to 1 Sv: <ul style="list-style-type: none"> photons <ul style="list-style-type: none"> - Hp(10) ±30 % from 15 keV to 9 MeV ²⁾ - Hp(0.07) ±30 % from 6 keV to 9 MeV beta <ul style="list-style-type: none"> - Hp(0.07) -20...+35% from 0.24 MeV to 0.80 MeV (E_{mean}) • angular response for photons: <ul style="list-style-type: none"> - Hp(10) ±20 % up to 60° at 65 keV - Hp(0.07) ±20 % up to 60° at 65 keV <p>¹⁾When calibrated after every 10 Sv of accumulated dose ²⁾ PTB approved model up to 7 MeV</p>
<p>Functional Characteristics</p>	<ul style="list-style-type: none"> • recording of official Hp(10) and Hp(0.07) doses • memory: <ul style="list-style-type: none"> - calibration date - dose reset dates - user identification
<p>Mechanical Characteristics</p>	<ul style="list-style-type: none"> • size: 41 x 44 x 12 mm, with holder 47 (95 with strap) x 49 x 13 mm (1.61 x 1.73 x 0.47 in , with holder 1.85 (3.74 with strap) x 1.93 x 0.51 in) • weight: 25 g (0.88 oz) , with holder 43 g (1.52 oz) • beta window: aluminized PI (app. 7 mg/cm²) • holder: anodized aluminium
<p>Environmental Characteristics</p>	<ul style="list-style-type: none"> • temperature range from - 10°C to +50 °C (14 °F to 122 °F) • enclosure class: IP 67 (waterproof)

The DBR-1 and DBR-2 Readers are designed to read DIS-1 dosimeters assembled in the DDH-2 Snap-in Dosimeter Holders. To obtain the most recent dose value, the user simply plugs the dosimeter into the reader and the dose values are displayed in a few seconds.



Typical energy response of personal dose equivalent Hp(10)



Typical energy response of personal dose equivalent Hp(0.07)



www.mirion.com
20996029_DIS1_EN_C

Mirion Technologies (MGPI) Inc
5000 Highlands Parkway
Suite 150
Smyrna Georgia 30082
USA
T +1.770.432.2744
F +1.770.432.9179

Mirion Technologies (MGPI) SA
BP 1
F-13113 Lamanon
France
T +33 (0) 4 90 59 59 59
F +33 (0) 4 90 59 55 18

Mirion Technologies (RADOS) Oy
P.O. Box 506
FIN-20101 Turku
Finland
T +358 2 4684 600
F +358 2 4684 601

Mirion Technologies (RADOS) GmbH
Ruhrstrasse 49
D-22761 Hamburg
Germany
T +49 40 85193 0
F +49 40 85193 256